

New Reservoir developed in pursuant of utmost usability.


PP PIPETTING RESERVOIR

20mL

Maximum volume
up to 50mL

NEW "Liquid Sweeper" has realized comfortable pipetting, developed focusing on usability and a feature includes the thought-through angle for gathering a liquid efficiently.

Liquid Sweeper




8 grooves, a very unique feature called "Liquid Sweeper", transfer a liquid drop from the top of the tip quickly.



Drop-on Cover



The most suitable angle for tip touch



The calculated slope angle gathers a liquid almost without its dead volume.

Convenient graduation



Easy-to-wash round corners



Autoclavable

Note: Minor deformation may occur due to thermal contraction, but it does not affect usage.

Wide size to cover 12ch pipette

Special cover for the reservoir to prevent dust or contamination. (Available only for sterilized type<1791-050-S>)

Note: The cover does not have sealability.

Cat. No.	Item	Unit	Case
1790-050-30	20mL, PP pipetting reservoir	30pcs.	10units
1791-050-S	20mL, PP pipetting reservoir, individual package (with cover)	Sterilized 10pcs.	5units