



204 short tip 5 plates/stack

10µL

705 standard tip 5 plates/stack

200µL



ECO PACK

Refill plates



Sterilized!

for ready to use.

- Environmentally-friendly
- Space-saving & Cost-saving
- Easy to refill



707 long tip 5 plates/stack

1000µL


10μL

Empty Rack
for ECO PACK 10 μL



1298-254RA

10 μL Short-tip ECO PACK



Real scale



Cat. No.	Item	Unit	Case
123R-254CS	10 μL, Short, Graduated, Natural, ECO PACK Sterilized	96tips×5plates×10stacks	5units
1298-254RA	Empty Rack for ECO PACK 10 μL	10racks	5units


200μL

Empty Rack
for ECO PACK 200 μL



1298-755RA

200 μL Standard-tip ECO PACK



Real scale



Cat. No.	Item	Unit	Case
123R-755CS	200 μL, Graduated, Natural, ECO PACK Sterilized	96tips×5plates×10stacks	5units
1298-755RA	Empty Rack for ECO PACK 200 μL	10racks	5units


1000μL

Empty Rack
for ECO PACK 1000 μL



1298-757RA

1000 μL Long-tip ECO PACK



Real scale



Cat. No.	Item	Unit	Case
123R-757CS	1000 μL, Long, Graduated, Natural, ECO PACK Sterilized	96tips×5plates×5stacks	5units
1298-757RA	Empty Rack for ECO PACK 1000 μL	5racks	—

Easy
to refill



Lift a plate



Place it into a rack



Refilling tips is completed



Nuclease-free certificate

This mark shows RNase/DNase/Human-DNA free. Products with this mark are manufactured in the cleanroom under our high quality control, and inspected to meet our strict quality standards. RNase/DNase-free certificates can be issued for each lot upon your request.

Something Different.
WATSON® BIO LAB
MADE IN JAPAN SINCE 1988

<https://watsonbiolab.com>

Sales agency:
WATSON CO., LTD.

Export Division:
2-2-7 Murotani Nishi-ku Kobe, 651-2241, JAPAN
TEL 81-78-991-4489 FAX 81-78-991-4491

Head Office:
1F, 14-17 Daikanyama-cho, Shibuya-ku, Tokyo,
150-0034, JAPAN
TEL 81-3-6615-3591 FAX 81-3-6427-0740

E-mail:
tcr@watson.co.jp

Manufacturer:
FUKAOKASEI CO., LTD.

※The specification and the price etc. of the product might change without a previous notice.