

New Reservoir developed in pursuant of utmost usability.

# PP PIPETTING RESERVOIR




**20mL**

Maximum volume  
up to 50mL

NEW "Liquid Sweeper" has realized comfortable pipetting, developed focusing on usability and a feature includes the thought-through angle for gathering a liquid efficiently.

## Liquid Sweeper




8 grooves, a very unique feature called "Liquid Sweeper", transfer a liquid drop from the top of the tip quickly.



## Drop-on Cover



### The most suitable angle for tip touch



The calculated slope angle gathers a liquid almost without its dead volume.

### Convenient graduation



### Easy-to-wash round corners



### Autoclavable

Note: Minor deformation may occur due to thermal contraction, but it does not affect usage.

### Wide size to cover 12ch pipette

Special cover for the reservoir to prevent dust or contamination. (Available only for sterilized type<1791-050-S>)

Note: The cover does not have sealability.

Cat. No.	Item	Unit	Case
1790-050-30	20mL, PP pipetting reservoir	30pcs.	10units
1791-050-S	20mL, PP pipetting reservoir, individual package (with cover)	<b>Sterilized</b> 10pcs.	5units