


Tube Racks

Freezer Storage Rack

- The storage box made of the polycarbonate. (Plate: Made of polypropylene)
- Excellent in heatproof and crash-proof.
- Conventional frozen preservation tube (to 2.0mL) can be applied to and the plate can be piled up.
- Number and alphabet for identification are marked to the partition.
- Possible to store in the liquid nitrogen (vapor phase) and in the deep freezing condition.

Cat.No.	Item	Dimensions(mm)	Packaging
1577-550	□ PC freezer storage rack, 50well	146x76x55	10pcs/Unit 10units/Case

□ Color : (W) white (B) blue (R) red (Y) yellow (G) green




Cat.No.	Item	Dimensions(mm)	Packaging
1577-850	□ PC freezer storage rack, 81well	131x131x52	10pcs/Unit 10units/Case

□ Color : (W) white (B) blue (R) red (Y) yellow (G) green

Piccolo Rack

- Mini tube rack for 24 tubes.
- Autoclavable rack made of the polycarbonate.
- Possible to store in the liquid nitrogen (vapor phase) and the deep freezing condition.



The bottom of the tube is clearly visible while setting up. Possible to put and store the crushed ice in this rack to cool the sample.





The notch(part of the furr) is made both on the body of the Watson Screw Cap Tube and on each of well of the Piccolo Rack so that the Screw Cap Tube is completely fixed. Thus, the cap of the Screw Tube can be easily opened and closed with one hand (much easier to open by using Screwcap Opener<#175-15C>).




Cat.No.	Item	Dimensions(mm)	Packaging
1522-2451	PC rack, 24well, for screwcap tubes	115x82x52	10pcs/Unit 5units/Case
1522-2453	PC rack, 24well, for micro tubes	115x82x52	10pcs/Unit 5units/Case

ECO Rack

- Ecology and economy.
- Made of 100% polypropylene recycled products.
- 96-well rack for 1.5/2.0mL, one size larger than 80 well rack.
- Writing space available on the side of the rack.




96well rack carrying two oblique 48well plates.



18mm pitch well interval. Multichannel pipette that installs the tip by skipping by one channel can be used.

Cat.No.	Item	Dimensions(mm)	Packaging
1521-902	ECO Rack, 96well	240.2x76.7x26.5	100pcs/Case


Color : (W) white (B) blue (Y) yellow (BL) black



The bottom of each well is shaped to stabilize the tube.

Microtube Preparation Rack

- Total 20 tubes of 1.5mL can be stored in a stand of two-row-stair shape.
- The bottom of tube is clearly visible while setting up.
- Anti-slip rubber equipped in two corners.



Cat.No.	Item	Dimensions(mm)	Packaging
1511-205	Microtube Preparation Rack	210x65x35	10pcs/Unit 5units/Case

PCR Tube Rack

- 96 of 0.2mL PCR single tubes can be stored. Frozen storage at -80°C. Autoclavable.
 - The black-colored shading type is recommended for the storage of the light sensitive sample.
 - Top and bottom is made of polycarbonate (black-colored shading type: made of polypropylene).
- Plate: Made of polypropylene



Cat.No.	Item	Dimensions(mm)	Packaging
1512-952	PCR Tube Rack (PC), 96well	118x82x40	10pcs/Unit 10units/Case
1512-952BL	PCR Tube Rack (PP), 96well, Black	118x82x40	10pcs/Unit 10units/Case

Something Different.
WATSON BIO LAB
MADE IN JAPAN SINCE 1988
<https://watsonbiolab.com>


Sales agency:
WATSON CO., LTD.

Export Division:
 2-2-7 Murotani Nishi-ku Kobe, 651-2241, JAPAN
 TEL 81-78-991-4489 FAX 81-78-991-4491

Head Office:
 1F, 14-17 Daikanyama-cho, Shibuya-ku, Tokyo,
 150-0034, JAPAN
 TEL 81-3-5615-3591 FAX 81-3-6427-0740

E-mail:
tcr@watson.co.jp

Manufacturer:
FUKAOKASEI CO., LTD.



*The specification and the price etc. of the product might change without a previous notice.