

Others

Disposable Loops / Needles

Pipetting Reservoir

Sterilized Scoop / Spoon



100mL

5mL

20mL

10mL


Disposable Loops / Needles

- Disposable inoculation loop and needle made of polystyrene.
- Comfortable and soft-touch handling.
- Available in 3 kinds of package: 1 piece, 10 pieces, 25 pieces.

Note: 10 pieces type and 25 pieces type are packed in a zipper bag.

Real scale

Material: Polystyrene



A convenient zipper bag used for 10 pcs / 25 pcs package.





Cat. No.	Item		Unit	Case
171-009-S	Disposable Inoculation Needle, Red	Sterilized	25 pieces × 40 bags	5 units
171-106-S1	Disposable Inoculation Loop, 1 µL, White, Individually wrapped package	Sterilized	250 pieces	5 units
171-109-S10	Disposable Inoculation Loop, 1 µL, White, 10 pcs. package	Sterilized	10 pieces × 100 bags	5 units
171-109-S	Disposable Inoculation Loop, 1 µL, White, 25 pcs. package	Sterilized	25 pieces × 40 bags	5 units
171-406-S1	Disposable Inoculation Loop, 10 µL, Blue, Individually wrapped package	Sterilized	250 pieces	5 units
171-409-S10	Disposable Inoculation Loop, 10 µL, Blue, 10 pcs. package	Sterilized	10 pieces × 100 bags	5 units
171-409-S	Disposable Inoculation Loop, 10 µL, Blue, 25 pcs. package	Sterilized	25 pieces × 40 bags	5 units



Pipetting Reservoir made of polypropylene 20 mL

User-Friendly Reservoir

Designed with an optimal slope for efficient liquid collection, ensuring comfortable pipetting. Autoclavable for convenience.



Liquid Sweeper

8 grooves, a very unique feature called "Liquid Sweeper", transfer a liquid drop from the top of the tip quickly.



Drop-on Lid

Special lid for the reservoir to prevent dust or contamination.

(Available only for sterilized type<1791-050-S>)

Note: The lid does not have sealability.


Dimensions: W117 mm × D51 mm × H6.2 mm

■ **The most suitable angle for tip touch**



The calculated slope angle gathers a liquid almost without its dead volume.

■ **Convenient graduation**



■ **Easy-to-wash round corners**




■ **Autoclavable**

Note: Minor deformation may occur due to thermal contraction, but it does not affect usage.

■ **Wide size to cover 12 ch pipette**

Pipetting Reservoir




Cat. No.	Item	Unit	Case
1790-050-30	Pipetting Reservoir, 20 mL, Bulk	30 pieces	10 units
1791-050-5	Pipetting Reservoir, 20 mL, with Lid, Individual package Sterilized	10 pieces	5 units


Disposable Doubled Individually Packaged

Sterilized Scoop/Spoon

- Perfect for sampling powders or semi-solid materials
- Double individual packaging is optimal for Clean Room use
- Inside scale for easy and precise sample collection
- Self-standing flat design




RNase / DNase / Human-DNA free
Products with this mark are manufactured in the clean room under our high quality control, and inspected to meet our strict quality standards. RNase/DNase-free certificates can be issued for each lot upon your request.



Pyrogen test
Results : Endotoxin level less than 0.05 EU/mL
* according to Internal test

The advantages of doubled individually package (e.g., at clean room)



Pre-clean room
Taking outer package out at pre-clean room, then moving to the second room, pass box. Any contaminants or materials on the inner package will be removed in these two rooms.

Clean Room
Taking Scoop/Spoon out from the inner package just before use prevent further contaminants or materials.

Package Material (Outer / Inner): Polyethylene

Sterilized Scoop

sterilized SCOOP





782-1005

100 mL

Length Overall: 204.6 mm
Maximum width of the scoop portion: 59.7 mm
Height: 31 mm



Material: polypropylene
(A meet the requirements of FDA)



Cat. No.	Item		Unit
782-1005	100 mL, Sterilized Scoop	Sterilized	50 pieces (individual packaging)

Shape and specifications are subject to change without notice for improvement.

Sterilized Spoon

Easy-insert-designed Slim head

sterilized SPOON



1 mL

782-001S
Length Overall: 157 mm
Maximum width of the spoon portion: 16.9 mm
Height: 9 mm



5 mL

782-005S
Length Overall: 184.9 mm
Maximum width of the spoon portion: 24.8 mm
Height: 9 mm




10 mL

782-010S
Length Overall: 185 mm
Maximum width of the spoon portion: 29.8 mm
Height: 18 mm




20 mL

782-020S
Length Overall: 207.7 mm
Maximum width of the spoon portion: 34.7 mm
Height: 18 mm



Measuring scale marked inside the spoon



Self-standing and Flat-bottom design



Designed to be ← horizontal when placed



Material: polypropylene
(A meet the requirements of FDA)



Cat. No.	Item		Unit
782-001S	1 mL, Sterilized Spoon	Sterilized	100 pieces (individual packaging)
782-005S	5 mL, Sterilized Spoon	Sterilized	100 pieces (individual packaging)
782-010S	10 mL, Sterilized Spoon	Sterilized	100 pieces (individual packaging)
782-020S	20 mL, Sterilized Spoon	Sterilized	100 pieces (individual packaging)

Shape and specifications are subject to change without notice for improvement.

Sterilized Scoop / Sterilized Spoon

